

Phonak Paradise

Features at a Glance



A Sonova brand

PHONAK
life is on

Phonak Paradise Models

Audéo™ P

Audéo P-R

P90/P70/P50
Rechargeable
Lithium-ion battery
Direct connectivity to most bluetooth enabled devices
Tap Control



Audéo P-RT

P90/P70/P50
Rechargeable, with telecoil
Lithium-ion battery
Direct connectivity to most bluetooth enabled devices
Tap Control



Audéo P-312

P90/P70/P50
Direct connectivity to most bluetooth enabled devices



Audéo P-13T

P90/P70/P50
Telecoil
Direct connectivity to most bluetooth enabled devices



Wireless Communication Portfolio Paradise Platform



1 TV Connector

Connect easily to a TV or stereo system. Turns Paradise hearing aids into wireless headphones with excellent stereo sound quality



2 Phonak Remote Control

Easily adjust volume and programs with the new physical remote control



3 Phonak PartnerMic™

A directional microphone for lapel use, allowing the hearing aid user to better communicate one-to-one in moderate noise

Rechargeable Charging Options



Phonak Mini Charger Case

- A compact charger for two Audéo M-R Audéo M-RT hearing aids
- Standard USB-C charging port
- Compatible with Audéo P-R & Audéo P-RT



Phonak Charger Case Combi

- A charger, drying kit and protective hard case all-in-one
- Includes a cleaning tool in the robust, splash-proof housing
- Compatible with Audéo P-R & P-RT



Phonak Power Pack

- For charging rechargeable hearing aids on the go
- Compatible with the Phonak Charger Case Combi only

Paradise Housing Colours

Hair & Skin



INSTANT FIT

Domes




CUSTOM EAR PIECES



Acoustic Coupling

S Receiver


Gain 46(2cc) • MPO 111(2cc)



Mild to moderate hearing loss

M Receiver


Gain 50(2cc) • MPO 114(2cc)



Mild to moderately-severe hearing loss

P Receiver


Gain 58(2cc) • MPO 122(2cc)



Moderate to severe hearing loss

UP Receiver

Gain 67(2cc) • MPO 122(2cc)



Moderate to severe to profound hearing loss

Discover the Wonder of Sound

There's nothing like the sound of Paradise



Unrivaled sound quality

- New sound processing chip - PRISM
- New motion sensor
- New APD 2.0
- New AutoSense OS 4.0

Personalized digital solutions

- New myPhonak 4.0 app
- New Hearing Screener 2.0

Universal connectivity

- Connects to smartphones, TV and more
- Made for all
- Multiple Bluetooth connections
- New Tap Control for easy access to Bluetooth functionalities

Premium P90

	FEATURES	BENEFITS	
AutoSense OS 4.0	 Calm Situation	Provides listening comfort and clarity in quiet, one on one conversation.	
	 Speech in Noise	Adapts and reduces noise from behind in real time and applies and optimizes the signal-to-noise ratio (SNR) in noisy environments.	
	 Speech in Loud Noise	Automatically helps you focus on a single speaker directly in front of you in an extremely noisy crowd, for eg. a party or restaurant.	
	 Speech in Car	Reduces car and road noise to make listening in the car more enjoyable.	
	 Comfort in Noise	For listening comfort in loud situations when no speech is present.	
	 Comfort in Echo	Reduces disturbing echo while letting you hear and understand speech naturally.	
	 Music	Offers a full and rich listening experience so you can enjoy your favourite music again.	
Connectivity	 Direct Connectivity	Automatic streaming to both ears from any bluetooth enabled device (phone calls, music, videos, e-books, podcasts and much more). Distinguishes between streamed speech and music signals, and automatically adjusts for an optimized sound quality.	
Manual	 Speech in 360°	Once activated, this automatically tracks and focuses on the dominant speaker, when you cannot face them, for e.g. driving with passengers behind or beside you.	
	 EchoBlock	Identifies signals blurred by echo/reverberations and removes them, providing you with clearer hearing.	
	 SoundRelax	Cushions sudden loud or unpleasant impulse sounds such as clanging dishes.	
	 WindBlock	Suppresses wind noise so you can enjoy spending time outdoors.	
	 WhistleBlock 20 channels	Elimination of annoying and embarrassing whistling from the hearing aids before it becomes audible to you or people around you.	
	 NoiseBlock 20 channels	Automatically reduces annoying background noise like engines or fans without reducing the voices you want to hear.	
	 QuickSync	Control both your instruments (program and/or volume changes) with one simple touch.	
Blocks	 UltraZoom Premium 33 channel adaptive directional mic	In noisy situations this feature zooms in on voices that are facing you. Noise from the side and back is reduced, for e.g. home activities or in meetings.	
	 Speech in Loud Noise	Automatically helps you focus on a single speaker directly in front of you in an extremely noisy crowd, for eg. a party or restaurant.	
	 SoundRecover2	Makes high pitched sounds easier to hear, for e.g. women's voices.	
	 Tinnitus Balance	Every Paradise model includes a "built-in" feature designed to reduce the annoyance of your tinnitus, especially in quiet environments.	
	 Dynamic Noise Cancellation	Automatically reduces background noise around you so that it is easier to hear the speech in front of you.	
	 Speech Enhancer	Designed to boost soft speech and speech at a distance in quiet situations.	
	 Tap Control (Audéo P-R & Audéo P-RT only)	Gives you easy access to voice assistant applications. This feature also allows you to answer and end phone calls from your smartphone and pause and resume streaming from your connected devices.	
	 Motion Sensor Hearing (Audéo P-R & Audéo P-RT only)	Detects whether you are standing or moving and adjusts the microphones to optimize your ability to hear speech.	
	Features		



Advanced P70

	FEATURES	BENEFITS
AutoSense OS 4.0	 Calm Situation	Provides listening comfort and clarity in quiet, one on one conversation.
	 Speech in Noise	Adapts and reduces noise from behind in real time and applies and optimizes the signal-to-noise ratio (SNR) in noisy environments.
	 Comfort in Noise	For listening comfort in loud situations when no speech is present.
	 Music	Offers a full and rich listening experience so you can enjoy your favourite music again.
Connectivity	 Direct Connectivity	Automatic streaming to both ears from any bluetooth enabled device (phone calls, music, videos, e-books, podcasts and much more). Distinguishes between streamed speech and music signals, and automatically adjusts for an optimized sound quality.
Manual	 Speech in Loud Noise	Helps you focus on a single speaker directly in front of you in an extremely noisy crowd, for eg. a party or restaurant.
	 Speech in 360°	Once activated, this automatically tracks and focuses on the dominant speaker, when you cannot face them, for e.g. driving with passengers behind or beside you.
Blocks	 SoundRelax	Cushions sudden loud or unpleasant impulse sounds such as clanging dishes.
	 WindBlock	Suppresses wind noise so you can enjoy spending time outdoors.
	 WhistleBlock 16 channels	Elimination of annoying and embarrassing whistling from the hearing aids before it becomes audible to you or people around you.
	 NoiseBlock 16 channels	Automatically reduces annoying background noise like engines or fans without reducing the voices you want to hear.
	 QuickSync	Control both your instruments (program and/or volume changes) with one simple touch.
Features	 UltraZoom Advanced 12 channel adaptive directional mic	Automatically zooms in on voices in noisy situations. Noise from the side and back is reduced, for eg. home activities or in meetings.
	 Speech in Loud Noise	Automatically helps you focus on a single speaker directly in front of you in an extremely noisy crowd, for eg. a party or restaurant.
	 SoundRecover2	Makes high pitched sounds easier to hear, for e.g. women's voices.
	 Tinnitus Balance	Every Paradise model includes a "built-in" feature designed to reduce the annoyance of your tinnitus, especially in quiet environments.
	 Tap Control (Audéo P-R Et Audéo P-RT only)	Gives you easy access to voice assistant applications. This feature also allows you to answer and end phone calls from your smartphone and pause and resume streaming from your connected devices.
	 Motion Sensor Hearing (Audéo P-R Et Audéo P-RT only)	Detects whether you are standing or moving and adjusts the microphones to optimize your ability to hear speech.

Standard P50

	FEATURES	BENEFITS
AutoSense OS 4.0	 Calm Situation	Provides listening comfort and clarity in quiet, one on one conversation.
	 Speech in Noise	Adapts and reduces noise from behind in real time and applies and optimizes the signal-to-noise ratio (SNR) in noisy environments.
	 Comfort in Noise	For listening comfort in loud situations when no speech is present.
Connectivity	 Direct Connectivity	Automatic streaming to both ears from any bluetooth enabled device (phone calls, music, videos, e-books, podcasts and much more). Distinguishes between streamed speech and music signals, and automatically adjusts for an optimized sound quality.
Manual	 Music	Offers a full and rich listening experience so you can enjoy your favourite music again.
Blocks	 WhistleBlock 12 channels	Elimination of annoying and embarrassing whistling from the hearing aids before it becomes audible to you or people around you.
	 NoiseBlock 12 channels	Automatically reduces annoying background noise like engines or fans without reducing the voices you want to hear.
Features	 QuickSync	Control both your instruments (program and/or volume changes) with one simple touch.
	 UltraZoom Standard Single channel adaptive directional mic	In noisy situations this feature zooms in on voices that are facing you. Noise from the side and back is reduced, for e.g. home activities or in meetings.
	 SoundRecover2	Makes high pitched sounds easier to hear, for e.g. women's voices.
	 Tinnitus Balance	Every Paradise model includes a "built-in" feature designed to reduce the annoyance of your tinnitus, especially in quiet environments.
	 Motion Sensor Hearing (Audéo P-R Et Audéo P-RT only)	Detects whether you are standing or moving and adjusts the microphones to optimize your ability to hear speech.

Life is on

At Phonak, we believe that well-hearing equates to well-being and is essential to living life to the fullest. For more than 70 years, we have remained passionate about creating a world where 'life is on' for everyone. Our innovative hearing solutions are designed for people of all ages and all degrees of hearing loss, to connect socially, thrive mentally and emotionally.

www.phonak.ca



Roger Installer

Roger installer is used alongside the Roger X receiver to install RogerDirect™ into any Paradise Hearing aid



Roger X Receiver

Required for RogerDirect™, the Roger X is the miniature universal Roger Receiver which is compatible with all Roger microphones and can be wirelessly installed into any Paradise hearing aid*



Roger Select

Roger Select is a versatile microphone ideal for situations where background noise is present such as restaurants



Roger TableMic II

Roger TableMic II is a microphone dedicated for working adults who participate in various meetings



Roger Pen

Bluetooth Connectivity providing an all-inclusive listening experience



Roger Easy Pen

Full Roger performance without Bluetooth functionality or manual microphone control

*Only Roger X with serial number with first four digits 1745-xxxx and higher

The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. Apple, the Apple logo, iPhone, and iOS are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google Inc.