

HEARING SHOULDN'T BE HARD WORK

RELY ON REXTON FOR A SMOOTH
ALL-ROUND HEARING EXPERIENCE



REXTON

OUR LATEST INNOVATIONS FOR RELIABLE HEARING



BiCore technology




BiCore hearing aids



MotionCore technology



M-Core hearing aids



Hearing aid model overview

BICORE B-Li M Rugged 8

Durability and a smooth listening experience



BICORE R-Li, R-Li T and R312 14

Our best Rexton RIC yet, now with BiCore technology



BICORE SR 18

Our smallest and lightest Slim RIC yet, now with BiCore technology



BICORE Custom-Li 20

Our first rechargeable Custom with the latest BiCore technology



M-CORE B-Li M, P, HP 24

Rechargeable BTE power, with MotionCore technology



M-CORE iX-CIC 31

Our smallest instant-fit CIC, with MotionCore technology



Portfolio overview 32

Proven technology

BRAND PROMISE

Our brand makes a solid statement – a reassurance to people with hearing loss that they can Rely on Rexton. No matter what life throws at them, we won't let them down.

RELY ON REXTON

CORE IDEA

Lifeproof is the idea that is the essence of the Rexton brand and it is present in everything we do. It is the takeaway we want people to remember and the most relevant reason for them to buy our products. It is a simple yet powerful message that differentiates us from the competition.

RELY ON
LIFEPROOF
REXTON



BRAND ATTRIBUTES



Proven
technology



Durable
construction



Practical and
easy-to-use



Dependable
company

CORE STORY

Rexton hearing aids get the job done every single time. They will perform reliably at home, at work, during exercise, and you can wear them in all weather. That's because they are built tough, using the highest quality materials, and tested rigorously to make sure they can answer all the challenges that life throws at you. Noisy rooms? No problem,

you can focus on what's being said. Drop them accidentally? They are impact proof. Caught in the rain? They are water resistant. Working in a harsh environment? Don't worry about it. And they can handle much, much more. Simply put, they'll keep you prepared for every situation. They are made to be Lifeproof.

BRAND MISSION

To deliver proven hearing technology in the most reliable way.

OUR AUDIENCE

Current and potential hearing aid users, and their families and Hearing Care Professionals (HCP).

BICORE TECHNOLOGY



Catching the important word of a colleague, registering the sound of the machinery, the people or the environment in any situation, is at the center of the Rexton BiCore technology.

Rexton is reliable and solid listening, even in demanding situations thanks to the breakthrough BiCore Technology – getting the job done is what is most important.

There are three special reasons why you can rely on every hearing aid in the range:

CLEAN SOUND

BiCore uses innovative sound processing to faithfully reconstruct the environmental soundscape for the hearing aid wearer.

The robust input stage and new A/D converter reliably handle sounds of all levels without distortion, ensuring that what comes in also goes out the same way, clean and smooth.

RELIABLE TECHNOLOGY

Patented Speech Preservation Technology® (SPT) minimizes unwanted noise and preserves speech sound.

SPT takes advantage of the spatial separation between target speech sounds and competing noise. This directional analysis allows for separate processing of the speech and noise before they are mixed together.

An ultra-quick wireless binaural link ensures left and right hearing aids are always synchronized and react instantly and precisely to changes in the sound environment.

LIFEPROOF

All BiCore hearing aids are Lifeproof with robust construction and proven design.

Rexton BiCore hearing aids are tested extensively to ensure they can withstand whatever life throws at them. They are durable enough to withstand everyday challenges and harsh environments.



MADE FOR IPHONE AND ANDROID (ASHA)

BiCore hearing aids can stream calls and audio from iPhones and Android phones.*



STREAMING

You may direct stream from your TV and other audio devices using the Smart Transmitter 2.4 or the Smart Mic.



TELECARE

Full remote fitting sessions via Connexx resulting in higher user convenience and satisfaction.



HEADSET MODE TAP CONTROL**

Rexton BiCore microphones get the job done so iOS users can enjoy hands free call.



MY VOICE 2.0**

General sound and own voice processed in parallel for a greater sound impression and user satisfaction.



AUTOECHO REDUCER**

Automated environment scan to reduce effects of disturbing reverberation on the perceived sound signal.

* Phones with Bluetooth 5.0+ and Android 10 (or higher) are capable of providing wireless streaming via ASHA, but ASHA must be enabled by each phone manufacturer and service carrier.

**Check performance level compatibility



BICORE B-LI M RUGGED

BICORE B-LI M RUGGED IS OUR TOUGHEST HEARING AID YET AND CAN HANDLE ALL THE KNOCKS AND SPILLS THAT CAN OCCUR IN EVERYDAY LIFE.

Based on BiCore technology for a smooth listening experience.




It has proven BiCore hearing technology that delivers a range of special features:

- > Ultra-fast adaption to changing environments
- > Patented speech preservation technology
- > Advanced binaural signal processing to reduce noise efficiently
- > Stream audio directly from your iOS or Android smartphone
- > Stream directly from your TV and other audio devices using the Smart Transmitter 2.4 or the Smart Mic

Unique tough building offering:

- > Water resistant construction
- > Sweat proof
- > Drop proof
- > Scratch resistant
- > ATEX certified

Colors



2 ccm Coupler*
 FOG (peak) 60 dB
 OSPL 90 (peak) 133 dB SPL

* All fitting range is based on 2 ccm coupler

TOUGH DESIGN IN A SUPER-COMPACT SIZE

To optimize toughness, a long list of enhancements have been made to the hearing aid structure. These include making the receiver airtight, sealing charging contacts, reinforcing receiver suspension, thickening the housing, and sealing all its parting lines. Rugged is also far smaller than conventional behind-the-ear hearing aids. Its dimensions are more comparable to those of a RIC.



Features / Performance Level	80	60	40
Channels/Controls/Programs	48/20/6	32/16/6	24/12/6
Soundpro 2.0	HighRes	HighRes	HighRes
My Voice 2.0	•	•	•
Direct Streaming	Made for iPhone / Android (ASHA)	Made for iPhone / Android (ASHA)	Made for iPhone / Android (ASHA)
Auto Volume	•	•	•
Wireless sync	•	•	•
Directionality	Automatic Adaptive, iOmni, Front/Back/Left/Right Automatic & Manual, Narrow	Automatic Adaptive, iOmni, Front/Back Automatic & Manual, Left/Right Manual, Narrow	Automatic Adaptive, iOmni, Narrow
Noise Reduction	Noise Management Sound Smoothing Directional	Noise Management Sound Smoothing Directional	Noise Management Sound Smoothing
Wind Noise Reduction	•	•	•
Reverb Reducer	•	•	-
AutoEcho Reducer	•	-	-
Bandwidth: Extension/Compression	•/•	-/•	-/•
Music Enhancer (presets)	3	3	1
Tinnitus	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy
Xphone	•	•	•
Acclimatization/Data Logging	•/•	•/•	•/•
Small earhook	○	○	○

Accessories & Apps

Charging Station B-M	Mandatory	Mandatory	Mandatory
Smart Key	○	○	○
Smart Transmitter 2.4	○	○	○
Smart Mic	○	○	○
Rexton App	○	○	○
BiCore CROS R312	○	○	○
BiCore CROS R-Li	○	○	○

Battery

Battery runtime (without streaming)	up to 39 h
Battery runtime (incl. 5 h streaming)	up to 36 h

• available - not available ○ optional

BICORE R-LI AND R-LI T

OUR BEST REXTON RIC YET, NOW WITH BICORE TECHNOLOGY, GIVES YOU DURABILITY AND A SMOOTH LISTENING EXPERIENCE.

The BiCore R-Li and R-Li T are rechargeable, they have Lifeproof microphones, 60% lower binaural latency rate (vs. the previous generation) and the lowest power consumption of any Rexton hearing aids.




- > Integrated motion sensor
- > Lithium-ion battery
- > Programmable rocker switch
- > Color Conversion Kits
- > Smart Demo device available
- > Programming via:
Noahlink Wireless

Colors



S-Receiver




Eartip 3.0 Open
 + Sleeve 3.0 Power
 + Earmold 3.0

2 ccm Coupler*

FOG (peak) 46 dB
OSPL 90 (peak) 110 dB SPL


M-Receiver



Eartip 3.0 Open
 + Sleeve 3.0 Power
 + Earmold 3.0

60 dB
119 dB SPL


P-Receiver



Sleeve 3.0 Power
 + Earmold 3.0

65 dB
122 dB SPL

HP-Receiver



Custom Shell (no vent)

75 dB
131 dB SPL

* All fitting range is based on 2 ccm coupler

Features / Performance Level	80	60	40	30	20
Channels/Controls/Programs	48/20/6	32/16/6	24/12/6	16/8/4	16/8/4
Soundpro 2.0	HighRes	HighRes	HighRes	HighRes	HighRes
My Voice 2.0	•	•	•	-	-
Direct Streaming	Made for iPhone / Android (ASHA)	Made for iPhone / Android (ASHA)	Made for iPhone / Android (ASHA)	Made for iPhone / Android (ASHA)	Made for iPhone / Android (ASHA)
Auto Volume	•	•	•	•	•
Wireless sync	•	•	•	•	•
Directionality	Automatic Adaptive, iOmni, Front/Back/Left/Right automatic & manual Narrow	Automatic Adaptive, iOmni, Front/Back/Left/Right manual, automatic & manual Narrow	Automatic Adaptive, iOmni, Narrow	Automatic Adaptive, iOmni	Automatic Adaptive, iOmni
Noise Reduction	Noise Management Sound Smoothing Directional	Noise Management Sound Smoothing Directional	Noise Management Sound Smoothing	Noise Management Sound Smoothing	Noise Management
Wind Noise Reduction	•	•	•	•	-
Reverb Reducer	•	•	-	-	-
AutoEcho Reducer	•	-	-	-	-
Bandwidth: Extension/Compression	• / •	- / •	- / •	- / •	- / •
Music Enhancer (presets)	3	3	1	1	-
Tinnitus	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	-
Xphone	•	•	•	•	•
Acclimatization/Data Logging	• / •	• / •	• / •	• / •	• / •
T-Coil	• (only R-Li T)	• (only R-Li T)	• (only R-Li T)	• (only R-Li T)	• (only R-Li T)
Accessories & Apps					
Smart Key	○	○	○	○	○
Smart Transmitter 2.4	○	○	○	○	○
Smart Mic	○	○	○	○	○
Rexton App	○	○	○	○	○
Travel Charger RIC/ Charging Station R/ Charging+ Station R	Mandatory	Mandatory	Mandatory	Mandatory	Mandatory
BiCore CROS R-Li	○	○	○	○	○
BiCore CROS R312	○	○	○	○	○
BiCore CROS SR	-	-	-	-	-
Battery**					
Battery runtime (without streaming)	up to 39 h (R-Li) up to 39 h (R-Li T)				
Battery runtime (incl. 5 h streaming)	up to 36 h (R-Li) up to 36 h (R-Li T)				

• available - not available ○ optional

** Battery runtime value based on S receiver.

BICORE R312

OUR BEST REXTON RIC YET, NOW WITH BICORE TECHNOLOGY, GIVES YOU DURABILITY AND A SMOOTH LISTENING EXPERIENCE.

Our newest battery RIC delivers sound that's true and clean. All have Lifeproof microphones, a 60% lower binaural latency rate (vs. the previous generation) and user-friendly rocker switch bumps for better usability and tactile feedback.




- > Integrated motion sensor
- > Battery size 312
- > Programmable rocker switch
- > Color Conversion Kits
- > Smart Demo device available
- > Programming via:
Noahlink Wireless

Colors



S-Receiver




Eartip 3.0 Open
 + Sleeve 3.0 Power
 + + Earmold 3.0

2 ccm Coupler*

FOG (peak) 46 dB
OSPL 90 (peak) 109 dB SPL


M-Receiver



Eartip 3.0 Open
 + Sleeve 3.0 Power
 + Earmold 3.0

58 dB
117 dB SPL


P-Receiver



+ Sleeve 3.0 Power
 Earmold 3.0

63 dB
120 dB SPL

HP-Receiver



Custom Shell (no vent)

74 dB
128 dB SPL

* All fitting range is based on 2 ccm coupler

Features / Performance Level	80	60	40
Channels/Controls/Programs	48/20/6	32/16/6	24/12/6
Soundpro 2.0	HighRes	HighRes	HighRes
My Voice 2.0	•	•	•
Direct Streaming	Made for iPhone / Android (ASHA)	Made for iPhone / Android (ASHA)	Made for iPhone / Android (ASHA)
Auto Volume	•	•	•
Wireless sync	•	•	•
Directionality	Automatic Adaptive, iOmni, Front/Back/Left/Right automatic & manual Narrow	Automatic Adaptive, iOmni, Front/Back/Left/Right manual, automatic & manual, Left/Right manual, Narrow	Automatic Adaptive, iOmni, Narrow
Noise Reduction	Noise Management Sound Smoothing Directional	Noise Management Sound Smoothing Directional	Noise Management Sound Smoothing
Wind Noise Reduction	•	•	•
Reverb Reducer	•	•	-
AutoEcho Reducer	•	-	-
Bandwidth: Extension/Compression	• / •	- / •	- / •
Music Enhancer (presets)	3	3	1
Tinnitus	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy
Xphone	•	•	•
Acclimatization/Data Logging	• / •	• / •	• / •
T-Coil	-	-	-

Accessories & Apps

Smart Key	○	○	○
Smart Transmitter 2.4	○	○	○
Smart Mic	○	○	○
Rexton App	○	○	○
BiCore CROS R-Li	○	○	○
BiCore CROS R312	○	○	○
BiCore CROS SR	-	-	-

Battery**

Battery runtime (without streaming)	up to 89 h
Battery runtime (incl. 20 h streaming)	up to 72 h

• available - not available ○ optional

** Battery runtime value based on S receiver.

BICORE SLIM RIC

GET THE PERFECT BLEND OF WEARING COMFORT AND STYLE WITH OUR SMALLEST AND LIGHTEST SLIM RIC YET, NOW WITH BICORE TECHNOLOGY.

Our newest Slim RIC is smaller and lighter than our previous models, has Lifeproof microphones, a 10% longer run time. Together with the Travel Charger it can provide up to 4 days of autonomy, longer than any other generation.



- > Integrated motion sensor
- > Lithium-ion battery
- > Smart Demo device available
- > Programming via:
Noahlink Wireless

Colors

Black | Silver



Black | Black Gloss



Cosmic Blue | Rose Gold




Snow White | Rose Gold



Snow White | Snow White Gloss




S-Receiver




Eartip 3.0 Open
Sleeve 3.0 Power
Earmold 3.0

M-Receiver



Eartip 3.0 Open
Sleeve 3.0 Power
Earmold 3.0

P-Receiver



Sleeve 3.0 Power
Earmold 3.0

2 ccm Coupler*

FOG (peak) 46 dB
OSPL 90 (peak) 110 dB SPL

60 dB
119 dB SPL

65 dB
122 dB SPL

* All fitting range is based on 2 ccm coupler

Features / Performance Level	80	60	40
Channels/Controls/Programs	48/20/6	32/16/6	24/12/6
Soundpro 2.0	HighRes	HighRes	HighRes
My Voice 2.0	•	•	•
Direct Streaming	Made for iPhone / Android (ASHA)	Made for iPhone / Android (ASHA)	Made for iPhone / Android (ASHA)
Auto Volume	•	•	•
Wireless sync	•	•	•
Directionality	Automatic Adaptive, iOmni, Front/Back/Left/Right automatic & manual Narrow	Automatic Adaptive, iOmni, Front/Back/automatic & manual, Left/Right manual, Narrow	Automatic Adaptive, iOmni, Narrow
Noise Reduction	Noise Management Sound Smoothing Directional	Noise Management Sound Smoothing Directional	Noise Management Sound Smoothing
Wind Noise Reduction	•	•	•
Reverb Reducer	•	•	-
AutoEcho Reducer	•	-	-
Bandwidth: Extension/Compression	•/•	-/•	-/•
Music Enhancer (presets)	3	3	1
Tinnitus	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy
Xphone	•	•	•
Acclimatization/Data Logging	•/•	•/•	•/•
T-Coil	-	-	-

Accessories & Apps

Smart Key	○	○	○
Smart Transmitter 2.4	○	○	○
Smart Mic	○	○	○
Rexton App	○	○	○
Slim-RIC Travel Charger 2 Charging+ Station SR	Mandatory	Mandatory	Mandatory
BiCore CROS R-Li	-	-	-
BiCore CROS R312	-	-	-
BiCore CROS SR	○	○	○

Battery**

Battery runtime (without streaming)	up to 20 h
Battery runtime (incl. 20 h streaming)	up to 17 h

• available - not available ○ optional ** Battery runtime value based on S receiver.

BICORE CUSTOM LI

OUR FIRST CUSTOM LI WITH THE LATEST BICORE TECHNOLOGY GIVES OPTIMAL DISCRETION AND HEARING PERFORMANCE

You can rely on our very first rechargeable custom model to deliver sound that's true and clean with Bluetooth streaming, contactless charging, iOS and Android compatibility and a sleek design.




- > Integrated motion sensor
- > Lithium-ion battery
- > Programmable push button
- > Programming via:
Noahlink Wireless


Colors




ITC 113/50




118/55




124/65



ITE 118/55



124/65



* All fitting range is based on 2 ccm coupler

Features / Performance Level	80	60	40
Channels/Controls/Programs	48/20/6	32/16/6	24/12/6
Soundpro 2.0	HighRes	HighRes	HighRes
Direct Streaming	Made for iPhone / Android (ASHA)	Made for iPhone / Android (ASHA)	Made for iPhone / Android (ASHA)
Auto Volume	•	•	•
Wireless sync	•	•	•
Directionality	Automatic Adaptive, iOmni, Front/Back/Left/Right automatic & manual Narrow	Automatic Adaptive, iOmni, Front/Back/automatic & manual, Left/Right manual, Narrow	Automatic Adaptive, iOmni, Narrow
Noise Reduction	Noise Management Sound Smoothing Directional	Noise Management Sound Smoothing Directional	Noise Management Sound Smoothing
Wind Noise Reduction	•	•	•
Reverb Reducer	•	•	-
AutoEcho Reducer	•	-	-
Bandwidth: Extension/Compression	•/•	-/•	-/•
Music Enhancer (presets)	3	3	1
Tinnitus	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy
Xphone	•	•	•
Acclimatization/Data Logging	•/•	•/•	•/•
T-Coil	-	-	-

Accessories & Apps

Smart Key	◦	◦	◦
Smart Transmitter 2.4	◦	◦	◦
Smart Mic	◦	◦	◦
Rexton App	◦	◦	◦
Custom Charger Station	Mandatory	Mandatory	Mandatory
BiCore CROS R-Li	◦	◦	◦
BiCore CROS R312	◦	◦	◦
BiCore CROS SR	-	-	-

Battery

Battery runtime (without streaming)	up to 28 h
Battery runtime (incl. 5 h streaming)	up to 24 h

• available - not available ◦ optional

MOTIONCORE TECHNOLOGY



Most hearing aids analyze hearing situations solely on acoustic criteria, then classify them into a small range of predefined acoustic environments. But when people rely on their senses to keep them safe and sound, hearing clearly in every situation is a must – especially when they're on the job and moving through the day, where environments can change quickly. Our M-Core hearing aids use proven MotionCore technology that automatically recognizes and adapts to different hearing situations. This ensures users can hear speech clearly and get the most reliable impression of their environment.



UP TO 20 TIMES MORE HEARING SITUATIONS AUTOMATICALLY ADAPTED THAN BEFORE



HIGH RESOLUTION SOUNDPRO

Using a combination of inputs including real-time motion data, signal to noise ratio, and sound directionality, High Resolution Soundpro can continuously adjust the sound settings to provide the most reliable impression of the user's environment, automatically.*



MOTION SENSOR

MotionCore technology features an integrated real-time motion sensor that detects a user's movement, and uses the data to even more accurately determine their sound situation. With this, M-Core hearing aids can automatically detect up to 127 different sound situations.



MY VOICE

The hearing aid recognizes the user's specific voice pattern, and integrates it naturally into what they hear. In turn, this provides the user with a natural-sounding voice that leads to high acceptance rates. And for that the user will communicate more and be more socially active, greatly improving their quality of life



TELECARE

Full remote fitting sessions via Connexx resulting in higher member convenience and satisfaction.



STREAMING

You may direct stream from your TV and other audio devices using the Smart Transmitter 2.4 or the Smart Mic.



REXTON ASSIST

Rexton Assist works on a live, learning network that continuously evolves to provide solutions to any challenge faced by the hearing aid wearer. Rexton Assist brings convenience for both the user and the dispenser, improving the hearing aid performance based on the needs of the wearer, at any time, in any place.


M-CORE B-LI M



VERSATILE AND RECHARGEABLE

Our popular BTE device features MotionCore technology and a rechargeable lithium-ion battery, for an accurate impression of your environment in every situation. Made for a wide range of hearing losses, M-Core B-Li M is designed and is available in Performance Levels 20-80. Your customers can choose from 10 colors, and have their device fitted open or closed depending on their needs.

- > Integrated motion sensor
- > Lithium-ion battery
- > Programmable rocker switch
- > ThinTube 3.0
- > Color Conversion Kits
- > Smart Demo device available
- > Programming via:
 - Noahlink Wireless
 - Rexton Assist



Colors



Features / Performance Level	80	60	40	30	20
Channels/Controls/Programs	48/20/6	32/16/6	24/12/6	16/8/4	16/8/4
Soundpro	HighRes	HighRes	HighRes	HighRes	HighRes
My Voice	•	•	•	-	-
Direct Streaming	Made for iPhone Smart Transmitter 2.4 and Smart Mic	Made for iPhone Smart Transmitter 2.4 and Smart Mic	Made for iPhone Smart Transmitter 2.4 and Smart Mic	Made for iPhone Smart Transmitter 2.4 and Smart Mic	Made for iPhone Smart Transmitter 2.4 and Smart Mic
Auto Volume	•	•	•	•	•
Wireless sync	•	•	•	•	•
Directionality	Automatic/Adaptive iOmni Front/Back Left/Right Narrow	Automatic/Adaptive iOmni Front/Back	Automatic/Adaptive iOmni	Automatic/Adaptive	Automatic/Adaptive
Noise Reduction	Noise Management Sound Smoothing Directional	Noise Management Sound Smoothing Directional	Noise Management Sound Smoothing	Noise Management Sound Smoothing	Noise Management
Wind Noise Reduction	Binaural	Binaural	Standard	Standard	-
Reverb Reducer	•	-	-	-	-
Bandwidth: Extension/Compression	• / •	- / •	- / •	- / •	- / •
Music Enhancer (presets)	3	3	1	1	-
Tinnitus	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	-
Xphone	•	•	•	-	-
Acclimatization/Data Logging	• / •	• / •	• / •	• / •	• / •
Small ear hook	○	○	○	○	○

Accessories & Apps	80	60	40	30	20
Charging Station	Mandatory	Mandatory	Mandatory	Mandatory	Mandatory
Smart Key	○	○	○	○	○
Smart Transmitter 2.4	○	○	○	○	○
Smart Mic	○	○	○	○	○
Rexton App	○	○	○	○	○
M-Core CROS R	○	○	○	-	-
M-Core CROS R-Li	○	○	○	-	-
M-Core CROS iX-CIC	-	-	-	-	-

Battery	80	60	40	30	20
Battery runtime (without streaming)			up to 24 h		
Battery runtime (incl. 5 h streaming)			up to 21 h		

• available - not available ○ optional


M-CORE B-LI P



POWERFUL AND RECHARGEABLE

M-Core B-Li P remains a powerful device that also offers the motion technology. Available in Performance Levels 20-80 for a wide range of applications.

- > Integrated motion sensor
- > Lithium-ion battery
- > Programmable rocker switch
- > ThinTube 3.0
- > Color Conversion Kits
- > Smart Demo device available
- > Programming via:
 - Noahlink Wireless
 - Rexton Assist



2 ccm coupler
 FOG (peak) 77 dB
 OSPL 90 (peak) 135 dB SPL

Colors

Black	Granite	Dark Champagne	Silver	Gray	Pearl White	Deep Brown	Sandy Brown	Rose Gold	Beige

Features / Performance Level	80	60	40	30	20
Channels/Controls/Programs	48/20/6	32/16/6	24/12/6	16/8/4	16/8/4
Soundpro	HighRes	HighRes	HighRes	HighRes	HighRes
My Voice	•	•	•	-	-
Direct Streaming	Made for iPhone Smart Transmitter 2.4 and Smart Mic	Made for iPhone Smart Transmitter 2.4 and Smart Mic	Made for iPhone Smart Transmitter 2.4 and Smart Mic	Made for iPhone Smart Transmitter 2.4 and Smart Mic	Made for iPhone Smart Transmitter 2.4 and Smart Mic
Auto Volume	•	•	•	•	•
Wireless sync	•	•	•	•	•
Directionality	Automatic / Adaptive iOmni Front & Back Left & Right Narrow	Automatic / Adaptive iOmni Front & Back	Automatic / Adaptive iOmni	Automatic / Adaptive	Automatic / Adaptive
Noise Reduction	Noise Management SoundSmoothing Directional	Noise Management SoundSmoothing Directional	Noise Management SoundSmoothing	Noise Management SoundSmoothing	Noise Management
Wind Noise Reduction	Binaural	Binaural	Standard	Standard	-
Reverb Reducer	•	-	-	-	-
Bandwidth: Extension/Compression	• / •	- / •	- / •	- / •	- / •
Music Enhancer (presets)	3	3	1	1	-
Tinnitus	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	-
Xphone	•	•	•	-	-
Acclimatization/Data Logging	• / •	• / •	• / •	• / •	• / •
T-Coil	•	•	•	•	•
Small ear hook	○	○	○	○	○

Accessories & Apps	80	60	40	30	20
Charging+ Station B-P / Charging Station B-P	Mandatory	Mandatory	Mandatory	Mandatory	Mandatory
Smart Key	○	○	○	○	○
Smart Transmitter 2.4	○	○	○	○	○
Smart Mic	○	○	○	○	○
Rexton App	○	○	○	○	○
M-Core CROS R	○	○	○	-	-
M-Core CROS R-Li	○	○	○	-	-
M-Core CROS iX-CIC	-	-	-	-	-

Battery	80	60	40	30	20
Battery runtime (without streaming)			up to 30 h		
Battery runtime (incl. 5 h streaming)			up to 27 h		

• available - not available ○ optional


M-CORE B-LI HP



EXTRA POWERFUL AND RECHARGEABLE

M-Core B-Li HP is a powerful device, suitable for those with severe to profound hearing loss, and combines the motion technology with rechargeability. Available in Performance Levels 20-80.

- > Integrated motion sensor
- > Lithium-ion battery
- > Programmable rocker switch
- > Color Conversion Kits
- > Smart Demo device available
- > Programming via:
 - Noahlink Wireless
 - Rexton Assist



2 ccm
coupler

FOG (peak)
82 dB

OSPL 90 (peak)
140 dB SPL

Colors



Features / Performance Level	80	60	40	30	20
Channels/Controls/Programs	34/20/6	32/16/6	24/12/6	16/8/4	16/8/4
Soundpro	HighRes	HighRes	HighRes	HighRes	HighRes
My Voice	•	•	•	-	-
Direct Streaming	Made for iPhone Smart Transmitter 2.4 and Smart Mic	Made for iPhone Smart Transmitter 2.4 and Smart Mic	Made for iPhone Smart Transmitter 2.4 and Smart Mic	Made for iPhone Smart Transmitter 2.4 and Smart Mic	Made for iPhone Smart Transmitter 2.4 and Smart Mic
Auto Volume	•	•	•	•	•
Wireless sync	•	•	•	•	•
Directionality	Automatic/Adaptive iOmni Front/Back Left/Right Narrow	Automatic/Adaptive iOmni Front/Back	Automatic/Adaptive iOmni	Automatic/Adaptive	Automatic/Adaptive
Noise Reduction	Noise Management Sound Smoothing Directional	Noise Management Sound Smoothing Directional	Noise Management Sound Smoothing	Noise Management Sound Smoothing	Noise Management
Wind Noise Reduction	Binaural	Binaural	Standard	Standard	-
Reverb Reducer	•	-	-	-	-
Bandwidth: Extension/Compression	•/•	-/•	-/•	-/•	-/•
Music Enhancer (presets)	3	3	1	1	-
Tinnitus	Sound Therapy	Sound Therapy	Sound Therapy	Sound Therapy	-
Xphone	•	•	•	-	-
Acclimatization/Data Logging	•/•	•/•	•/•	•/•	•/•
T-Coil	•	•	•	•	•
Small ear hook	○	○	○	○	○

Accessories & Apps	80	60	40	30	20
Charging+ Station B-HP / Charging Station B-HP	Mandatory	Mandatory	Mandatory	Mandatory	Mandatory
Smart Key	○	○	○	○	○
Smart Transmitter 2.4	○	○	○	○	○
Smart Mic	○	○	○	○	○
Rexton App	○	○	○	○	○
M-Core CROS R	○	○	○	-	-
M-Core CROS R-Li	○	○	○	-	-
M-Core CROS iX-CIC	-	-	-	-	-

Battery	80	60	40	30	20
Battery runtime (without streaming)				up to 61 h	
Battery runtime (incl. 5 h streaming)				up to 57 h	

• available - not available ○ optional




M-CORE SR

WITH EXCHANGEABLE RECEIVERS



Our discreet Slim-RIC design features proven MotionCore technology that automatically recognizes and adapts to up to 20 times more hearing situations than before.

- > Integrated motion sensor
- > Lithium-ion battery
- > Smart Demo device available
- > Programming via:
 - Noahlink Wireless
 - Rexton Assist



- Colors
- Black | Silver
 - Black | Black Gloss
 - Cosmic Blue | Rose Gold
 - Snow White | Rose Gold
 - Snow White | Snow White Gloss
- * Colors not available for all performance levels

Features / Performance Level	80	60	40	30	20
Channels/Controls/Programs	48/20/6	32/16/6	24/12/6	16/8/4	16/8/4
Soundpro	HighRes	HighRes	HighRes	HighRes	HighRes
My Voice	•	•	•	-	-
Direct Streaming	Made for iPhone Smart Transmitter 2.4 and Smart Mic	Made for iPhone Smart Transmitter 2.4 and Smart Mic	Made for iPhone Smart Transmitter 2.4 and Smart Mic	Made for iPhone Smart Transmitter 2.4 and Smart Mic	Made for iPhone Smart Transmitter 2.4 and Smart Mic
Auto Volume	•	•	•	•	•
Wireless sync	•	•	•	•	•
Directionality	Automatic / Adaptive iOmni Front & Back Left & Right Narrow	Automatic / Adaptive iOmni Front & Back	Automatic / Adaptive iOmni	Automatic / Adaptive	Automatic / Adaptive
Noise Reduction	Noise Management SoundSmoothing Directional	Noise Management SoundSmoothing Directional	Noise Management SoundSmoothing	Noise Management SoundSmoothing	Noise Management
Wind Noise Reduction	Binaural	Binaural	Standard	Standard	-
Reverb Reducer	•	-	-	-	-
Bandwidth: Extension/Compression	• / •	- / •	- / •	- / •	- / •
Music Enhancer (presets)	3	3	1	1	-
Tinnitus	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	-
Xphone	•	•	•	-	-
Acclimatization/Data Logging	• / •	• / •	• / •	• / •	• / •

Accessories & Apps	80	60	40	30	20
Slim-RIC Travel Charger	Mandatory	Mandatory	Mandatory	Mandatory	Mandatory
Smart Key	○	○	○	○	○
Smart Transmitter 2.4	○	○	○	○	○
Smart Mic	○	○	○	○	○
Rexton App	○	○	○	○	○
M-Core CROS R	-	-	-	-	-
M-Core CROS R-Li	-	-	-	-	-
M-Core CROS iX-CIC	-	-	-	-	-

Battery**	80	60	40	30	20
Battery runtime (without streaming)				up to 19 h	
Battery runtime (incl. 5 h streaming)				up to 15 h	

• available - not available ○ optional ** Battery runtime value based on S receiver.


M-CORE iX-CIC

WITH THE HIDDEN POWER OF MOTIONCORE



Our virtually invisible CIC design features proven MotionCore technology that automatically recognizes and adapts to up to 10 times more hearing situations than before. This ensures users hear speech clearer and get the most reliable impression of their environment, automatically. M-Core iX hearing aids are designed to be discreet, practical and easy-to-use.

- > 10 size battery
- > Smart Demo device available
- > Programming via:
ConnexxAir
Flex connector



Colors

Shell: Red Blue Mocha Black

Faceplate:

Features / Performance Level	80	60	40	30	20
Channels/Controls/Programs	48/20/6	32/16/6	24/12/6	16/8/4	16/8/4
Soundpro	HighRes	HighRes	HighRes	HighRes	HighRes
My Voice	-	-	-	-	-
Direct Streaming	-	-	-	-	-
Wireless sync	•	•	•	•	•
Directionality	Automatic/Adaptive iOmni Front	Automatic/Adaptive iOmni Front	Automatic/Adaptive iOmni	-	-
Noise Reduction	Noise Management SoundSmoothing Directional	Noise Management SoundSmoothing Directional	Noise Management SoundSmoothing	Noise Management SoundSmoothing	Noise Management
Wind Noise Reduction	Standard	Standard	Standard	Standard	-
Reverb Reducer	•	-	-	-	-
Bandwidth: Extension/Compression	• / •	- / •	- / •	- / •	- / •
Music Enhancer (presets)	3	3	1	1	-
Tinnitus	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	Sound Therapy Notch Therapy	-
Xphone	•	•	•	-	-
Acclimatization/Data Logging	• / •	• / •	• / •	• / •	• / •

Accessories & Apps	80	60	40	30	20
Smart Key	○	○	○	○	○
Rexton App	○	○	○	○	○
Smart Mic	-	-	-	-	-
Smart Transmitter 2.4	-	-	-	-	-
M-Core CROS iX-CIC	•	•	•	-	-

Battery	80	60	40	30	20
Battery runtime	up to 70 h				

• available - not available ○ optional

PORTFOLIO OVERVIEW

RELIABLE HEARING OPTIONS FOR YOUR CUSTOMERS

RIC

Slim RIC

BTE



Product

**BiCore
R-Li**

**BiCore
R-Li T**

**BiCore
R312**

BiCore SR

B-Li M Rugged

**M-Core
B-Li
M,P,HP**

Battery

Li-ion

Li-ion

312

Li-ion

Li-ion

Li-ion

Programming
Device

Noahlink Wireless

Noahlink Wireless

Noahlink Wireless

Noahlink Wireless

Noahlink Wireless

Noahlink Wireless

TeleCare

Connexx Full Live
Tuning

Connexx Full Live
Tuning

Connexx Full Live
Tuning

Connexx Full Live
Tuning

Connexx Full Live
Tuning

Connexx Full Live
Tuning

Bluetooth Capable

Yes

Yes

Yes

Yes

Yes

Yes

T-Coil
Option

No

Yes

No

No

Yes

Yes*

Rexton Assist

Yes

Yes

Yes

Yes

Yes

Yes

* not available for M products

Instant Fit

Custom



**M-Core
iX-CIC**

**BiCore
Custom Li**

10

Li-ion

Connexx Air, HiPro II

Noahlink Wireless

Portal
Basic Tuning

Connexx Full Live
Tuning

No

Yes

No

No

No

Yes



RELY ON REXTON

It's a tough world out there, and even tougher when your hearing starts to go. At Rexton, we won't let hearing loss hold you back. Good hearing is critical: From getting the job done, to getting home safely and being there for the people who count on you, we know what's at stake. That's why we work hard to deliver proven hearing technology in the most reliable way. We understand what reliability means in real life, and make practical, easy-to-use products you can count on. Our hearing aids have been getting the job done since 1955, so whatever life throws at you, you can rely on Rexton.

We believe that nothing hearing related should ever be left to chance, both for the people with hearing loss, and those who depend upon them. We understand their challenges and we meet them with proven hearing technology that is made Lifeproof so it can be counted upon to perform reliably in even the harshest environments. At work, at home, during leisure time, exercise and in all weather, you can always RELY ON REXTON.

**GETTING
THE JOB DONE
SINCE 1955**

[rexton.com](https://www.rexton.com)

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos GmbH is under license. Other trademarks and trade names are those of their respective owners. Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Google Play is a trademark of Google Inc. "Made for iPhone" and "Made for iPad" mean that an electronic accessory has been designed to connect specifically to iPhone or iPad, respectively, and has been certified by the developer to meet Apple performance standards. Apple is not responsible for the operation of this device or its compliance with safety and regulatory standards. Please note that the use of this accessory with iPhone or iPad may affect wireless performance. iPad, iPhone are trademarks of Apple Inc., registered in the U.S. and other countries.

Order No. 04538-99C1-7600 © 2022 Sivantos GmbH