

Signia app



### You're in control.

Enjoy your hearing aids to the full with just one app for everything you need:

- All controls are tailored to your personal preferences
- Audio streaming, remote control and remote support
- Discreetly control the volume and other settings while streaming audio directly into your ears
- Stay in touch with your Hearing Care Professional and get remote support when you need it



Compatible with Android and iOS. The app are free of charge.

Pure 312 X



**Life doesn't stand still. 90% of our activities involve motion. As we move, the acoustic environment changes and hearing can go from easy to nearly impossible in a matter of seconds.**

With Pure 312 X, you can move freely and still hear what's important in superb clarity. The world's first hearing aids with acoustic motion sensors enable you to hear clearly wherever you go and be an active part of every conversation. Pure 312 X delivers a hearing experience so natural, it's sure to move you.



The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice.

#### Legal Manufacturer

Signia GmbH  
Henri-Dunant-Strasse 100  
91058 Erlangen  
Germany

Order No. 03986-99C1-7600 | Printed in the UK | JN10020  
© 08.2019 Signia GmbH

Find us on:

[signia-hearing.co.uk](http://signia-hearing.co.uk)



[facebook.com/signiahearing](https://facebook.com/signiahearing)



[youtube.com/signiahearing](https://youtube.com/signiahearing)

Pure 312 X

**signia**

Life sounds brilliant.



The new Pure 312 X

**HEAR WHAT MATTERS TO YOU.**



[signia-hearing.co.uk/pure312x](http://signia-hearing.co.uk/pure312x)

**HEARING:  
IT'S A MOVING  
EXPERIENCE**



Pure 312 X

# SOUNDS JUST MADE FOR YOU



Everyone is different. And everyone has a unique hearing experience depending on where they are and what they are doing. But until now your hearing aids could only listen and make educated guesses about the type of situation you're in.



Pure 312 X redefines what modern hearing aids can do:

- Understand speech in every situation
- Hear your environment in greater detail than ever
- Enjoy crystal clear sound tailored to your activities
- Be independent with fully rechargeable and ultra-reliable hearing aids

**Hear precisely what matters to you.**

# ALL THE POWER YOU NEED

The new Pure 312 X delivers a powerful performance for an outstandingly personal hearing experience. All day long.



Pure 312 X

Pure 312 X with T-coil

## At a Glance

- Equipped with revolutionary acoustic-motion sensors
- Completely integrated T-coil
- Bluetooth connectivity for direct streaming of phone calls, music and TV
- With tinnitus treatment options

T-coil



## Enjoy the performance, with complete discretion.

For maximum enjoyment and convenience when you're out and about, Pure 312 X's T-coil function picks up the signal from compatible public sound systems, such as those found in theaters, and delivers it directly to your ears. Plus, it's completely integrated into the housing for an elegant and discreet solution. So you can enjoy the performance without standing out from the crowd.



## Stay connected.


In today's world, everything and everyone is connected. Your hearing aids are too. Thanks to their Bluetooth compatibility, Pure 312 X hearing aids let you stream sound from your favorite devices straight to your ears.

- Calls on your smartphone will be automatically streamed directly to your hearing aids
- Listen to your favorite music or TV programs via Bluetooth

## Ease tinnitus. The easy way.

For tinnitus sufferers, Pure 312 X features a range of therapy options to make your life more comfortable.

Available colours



Find out more about how Pure 312 X can give you a more personalised hearing experience at [signia-hearing.co.uk/pure312x](http://signia-hearing.co.uk/pure312x)

You will also find more information about useful accessories and the Signia app, or ask your hearing care professional.