

Welcome to the **Vivante**[™] *experience*

Five steps to a totally personalized hearing solution

Vivante Product Selection Guide 1.0



unitron[™] Love the experience

A Sonova brand

Five steps to the hearing solution that's right for you

Step one

Where do you experience hearing challenges?

Step two

Explore the technology levels

Step three

Pick a style that suits you

Step four


Color your choice

Step five

Connect to convenience

Step one

Where do you experience hearing challenges?



Step two

Explore the technology levels

Now we'll select the technology level that best matches your listening lifestyle.

What's included

What you're missing

Level 9 ★★★★★

Premium performance in 8 environments and our most sophisticated features that automatically let you focus on speech and know where it's coming from in crowds and loud background noise (where you spend ___% of your time). Enjoy the ability to hear in quiet as well as delightful music automatically. Plus engage in conversations in the car with a special program that pinpoints and amplifies speech from any direction no matter where you are sitting.

Nothing. You've got everything available to suit your unique lifestyle.

Level 7 ★★★★

Advanced performance in 6 environments with the ability to help you focus on speech coming from the front and sides in challenging situations, and hear in quiet situations, plus an automatic music program and a manual car program that can be pre-set to focus on speech from one location.

The most sophisticated features to focus on speech from different directions in noisy conversations in loud noise and the car.

Level 5 ★★★

Standard performance in 4 environments that includes features to help you in multiple quiet and noisy environments.

Automatic optimization for conversation in small groups and automatic engagement of a music program.

Level 3 ★★

The essentials in 2 environments. Limited ability to adjust automatically to accommodate more challenging listening situations. The flexibility to customize your experience with the fully featured Remote Plus app and control your streamed audio, calls and virtual assistants with tap control.

Automatic adjustment for comfort in quiet and noisy situations without speech.

Improve quiet and challenging conversations with our latest advanced technology.

Automatic program

8 environments



Conversation in quiet



Conversation in noise



Quiet



Noise



Conversation in a small group



Music



Conversation in a crowd



Conversation in loud noise

Manual programs



Conversation in loud noise



360 conversation in car

Various manual programs can be set

6 environments



Conversation in quiet



Conversation in noise



Quiet



Noise



Conversation in a small group



Music



Conversation in loud noise

Car with speech

Various manual programs can be set

4 environments



Conversation in quiet



Conversation in noise



Quiet



Noise

Various manual programs can be set

2 environments



Conversation in quiet



Conversation in noise


Various manual programs can be set

Sound performance by technology level

Low complexity

Level 9

Premium
★★★★★




Pinna Effect 2

Binaural microphone mode that is designed to compensate for cues typically lost with hearing aids, which are needed for localization.

Level 7

Advanced
★★★★




Pinna Effect 2

Binaural microphone mode that is designed to compensate for cues typically lost with hearing aids, which are needed for localization.

Level 5

Standard
★★★



Pinna Effect 2

Binaural microphone mode that is designed to compensate for cues typically lost with hearing aids, which are needed for localization.

Level 3

Essential
★★

Pinna Effect

Minimizes the impact of acoustic cues typically lost with hearing aids.

High complexity



Speech from front

AutoFocus 360 identifies and focuses on speech from the front to provide a greater increase of speech than other technology levels. Dynamic noise reduction is active and adjustable to further reduce surrounding noise. For face-to-face conversations, in loud noise, and when the wearer is stationary, HyperFocus engages within the conversation in loud noise environment.



Speech from side





Speech from back

AutoFocus 360 can identify and focus on speech from any direction while reducing interfering noise.



Speech from front

AutoFocus identifies, focuses on and increases speech from the front. Dynamic noise reduction is active and adjustable to further reduce surrounding noise.



Speech from side and back

AutoFocus applies a wide front beam to capture speech partially from the sides.



All high complexity situations

Applies a wide front beam capturing the speech from the front.



All high complexity situations

Applies a wide front beam capturing speech from the front.

Step three

Pick a style that suits you

Choose from products available
on the **Vivante™** platform.

Moxi™ Vivante rechargeable hearing aids are designed to focus on comfort, aesthetics and ease-of-use. Conveniently stream phone calls and other media, like music, podcasts and digital assistants right to your hearing aids, easily controlling them with just a tap.*



Moxi V-R

Our smallest rechargeable
with direct connectivity

Receiver-in-canal (RIC)



Moxi V-RT

Our rechargeable with direct
connectivity and telecoil

Receiver-in-canal (RIC)

Charging case for
Moxi V-R and Moxi V-RT
rechargeable
hearing aids



Products on the **Blu** platform

Other styles that could suit your hearing needs are available on the Blu platform**



Stride™ B-PR

Our rechargeable BTE with direct connectivity

Behind-the-ear (BTE)



Stride B-UP

Our Ultra Power BTE with direct connectivity, telecoil, and traditional 675 disposable battery

Behind-the-ear (BTE)

Charging case for
Stride rechargeable
hearing aids



Insera™ B-312

Our wireless ITE model with direct connectivity

In-the-ear (ITE)



Insera B-312 NW O

Our non-wireless omni ITE model with a telecoil

In-the-ear (ITE)



Insera B-10 NW O

Our smallest non-wireless omni ITE model with a telecoil

In-the-ear (ITE)

Hearing aids shown at actual size. Charging case not shown at actual size

*Tap control is available for rechargeable models only.

** Hearing aids on the Blu platform differ from the Vivante platform. Speak with your hearing care professional for more information.

Step four

Color your choice



Sand
(R8)



Amber
(P2)



Espresso
(P4)



Platinum
(P6)



Pewter
(P7)



Charcoal
(P8)



Beige
(26)



Tan
(14)



Cocoa
(22)



Brown
(28)



Black
(06)*

Colors shown on the Moxi V-R style are also available on all Moxi V-RT and Stride Blu styles.
Moxi V-R shown at actual size.

*Available for Inera B-312 and Inera B-10 IIC styles only.

Step five

Connect to convenience



TV Connector

Stream your favorite media wirelessly in stereo sound with the TV Connector.



PartnerMic

PartnerMic is an easy-to-use lapel microphone ideal for one-on-one conversations.



Remote Control

Discreetly adjust volume and switch programs with the Remote Control.

Remote Plus app

The Remote Plus app puts easy personalization at your fingertips.



Scan to download the Remote Plus app.



iOS



Android

Need more help in challenging listening situations like meetings or conversations in noise or over distance?
Let's discuss if a Roger™ device could be a good fit for you.

Accessories shown at actual size.

Accessories shown are available with all Vivante products and wireless products on the Blu platform.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

Google Play and the Google Play logo are trademarks of Google Inc

What's next?

- Discuss your needs

- Try at home

- Come in for a follow-up assessment

- Personalize your solution

- Come in for your final fitting on:

Recommendation:			
Technology level:	Style:	Warranty:	Cost:
Accessories:			

Notes:



At Unitron, we empower people with life-enhancing hearing experiences that fit seamlessly into their world. Our sound performance technology, Experience Innovations, and intuitive design work perfectly together for unmatched personalization and optimization. Because everyone deserves to **Love the experience.™**

© 2023 Unitron. All rights reserved.

028-6711-02/V1.00/2023-04

unitron.com



sonova
HEAR THE WORLD